



# Pristine Connect™: the different user types.

	Facility Admin	OR Staff	Physician Asst	Scheduler	Patient	Surgeon
Schedule a Case	✓	✓		✓		
Pair Summit to Patient	✓	✓				
View Images & Videos	✓	✓	✓		✓	✓
Share Images & Videos			✓			✓
Add Notes to Images & Videos			✓			✓
View Surgeon Case Schedule	✓		✓			✓
View Facility Case Schedule	✓	✓				
View Inventory Details	✓					
View Order Details	✓					
Add Surgeons or Staff	✓					
Edit Surgeons or Staff	✓					
Assign Staff Roles and Responsibilities	✓					
Mobile app						✓